

NPY (Human, Rat, Mouse)

| Catalog# | Standard size | Price | Order |
|----------|---------------|----------|--------------------------------------|
| 049-03 | 200 µg | \$105.00 | <input type="button" value="ORDER"/> |

Product description

More information

References

Normal Plasma Peptide/Protein Levels

Neuropeptide Y (NPY)

| | |
|----------------------------|---|
| Catalog # | 049-03 |
| Standard Size | 200 µg |
| Sequence | Tyr - Pro - Ser - Lys - Pro - Asp - Asn - Pro - Gly - Glu - Asp - Ala - Pro - Ala - Glu - Asp - Met - Ala - Arg - Tyr - Tyr - Ser - Ala - Leu - Arg - His - Tyr - Ile - Asn - Leu - Ile - Thr - Arg - Gln - Arg - Tyr - NH ₂ |
| Species | Human, Rat, Mouse |
| M.W | 4271.73 |
| Mass Spec Analysis | Exhibits correct M.W. |
| Solubility | Soluble in water. |
| Storage Condition | Store peptide in DRY form at 0-5°C. |
| Appearance | White powder |
| Content | The contents of this vial have been accurately determined to contain 200µg of NET peptide. |
| HPLC Analysis | Purity ≥96% |
| Additional articles | For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions. {Not Ovine.} |