

Catestatin (Human)

Catalog#	Standard size	Price	Order
053-27	200 µg	\$105.00	<input type="button" value="ORDER"/>

- Product description
- More information
- References

[Chromogranin A](#)
[Catestatin](#)
[Serpinin](#)

Catalog #	053-27
Standard Size	200 µg
Sequence	Ser - Ser - Met - Lys - Leu - Ser - Phe - Arg - Ala - Arg - Ala - Tyr - Gly - Phe - Arg - Gly - Pro - Gly - Pro - Gln - Leu
Species	Human
M.W	2326.71
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 200µg of NET peptide.
HPLC Analysis	Purity ≥90%