

Catestatin (Mouse)

Catalog#	Standard size	Price	Order
053-28	200 µg	\$105.00	<input type="button" value="ORDER"/>

- Product description
- More information
- References

[Chromogranin A](#)
[Catestatin](#)
[Serpinin](#)

Catalog #	053-28
Family	Please call for status of this item
Standard Size	200 µg
Sequence	Arg - Ser - Met - Arg - Leu - Ser - Phe - Arg - Thr - Arg - Gly - Tyr - Gly - Phe - Arg - Asp - Pro - Gly - Leu - Gln - Leu
Species	Mouse
M.W	2513.91
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 200µg of NET peptide.
HPLC Analysis	Purity ≥97%