

Neuronostatin-13 (Human, Porcine)

Catalog#	Standard size	Price	Order
060-50	200 µg	\$85.00	<input type="button" value="ORDER"/>

Product description

More information


References

Neuronostatin

Catalog #	060-50
Standard Size	200 µg
Sequence	Leu - Arg - Gln - Phe - Leu - Gln - Lys - Ser - Leu - Ala - Ala - Ala - Ala - NH ₂
Species	Human, Porcine
M.W	1415.67
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0 - 5°C. For best results and reproducibility, rehydrate peptide just before use. Do not re-freeze any unused portion.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 200 µg of NET peptide
HPLC Analysis	Purity ≥95%

Graphic Data

Mass Spec Analysis of Neuronostatin-13 (Human, Porcine)



HPLC Analysis of Neuronostatin-13 (Human, Porcine)

