

## Neuronostatin-19 (Human)

Catalog#	Standard size	Price	Order
060-53	100 µg	\$170.00	<input type="button" value="ORDER"/>

### Product description

### More information

### References

#### Neuronostatin

<b>Catalog #</b>	060-53
<b>Standard Size</b>	100 µg
<b>Sequence</b>	Ala - Pro - Ser - Asp - Pro - Arg - Leu - Arg - Gln - Phe - Leu - Gln - Lys - Ser - Leu - Ala - Ala - Ala - Ala - NH <sub>2</sub>
<b>Species</b>	Human
<b>M.W</b>	2039.37
<b>Mass Spec Analysis</b>	Exhibits correct M.W.
<b>Solubility</b>	Soluble in water.
<b>Storage Condition</b>	Store peptide in DRY form at 0 - 5°C.
<b>Appearance</b>	White powder
<b>Content</b>	The contents of this vial have been accurately determined to contain 100 µg of NET peptide
<b>HPLC Analysis</b>	Purity ≥95%
<b>Additional articles</b>	For best results and reproducibility, rehydrate peptide just before use. Do not re-freeze any unused portion.