

## Irisin (Human, Rat, Mouse, Canine)

Catalog#	Standard size	Price	Order
067-16	100 µg	\$450.00	<input type="button" value="ORDER"/>

### Product description

### More information

[Irisin - Antibody Stainings, Sequences & References](#)

[Irisin Synthesis and Its Applications](#)

<b>Catalog #</b>	067-16
<b>Standard Size</b>	100 µg
<b>Sequence</b>	Asp - Ser - Pro - Ser - Ala - Pro - Val - Asn - Val - Thr - Val - Arg - His - Leu - Lys - Ala - Asn - Ser - Ala - Val - Val - Ser - Trp - Asp - Val - Leu - Glu - Asp - Glu - Val - Val - Ile - Gly - Phe - Ala - Ile - Ser - Gln - Gln - Lys - Lys - Asp - Val - Arg - Met - Leu - Arg - Phe - Ile - Gln - Glu - Val - Asn - Thr - Thr - Thr - Arg - Ser - Cys - Ala - Leu - Trp - Asp - Leu - Glu - Glu - Asp - Thr - Glu - Tyr - Ile - Val - His - Val - Gln - Ala - Ile - Ser - Ile - Gln - Gly - Gln - Ser - Pro - Ala - Ser - Glu - Pro - Val - Leu - Phe - Lys - Thr - Pro - Arg - Glu - Ala - Glu - Lys - Met - Ala - Ser - Lys - Asn - Lys - Asp - Glu - Val - Thr - Met - Lys - Glu
<b>Species</b>	Human, Mouse, Rat, Canine
<b>M.W</b>	12587.25
<b>Mass Spec Analysis</b>	Exhibits correct M.W.
<b>Storage Condition</b>	Store peptide in DRY form at 0-5°C For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions
<b>Appearance</b>	White Powder
<b>Content</b>	100 µg NET peptide
<b>HPLC Analysis</b>	Purity ≥ 95%

Thursday 20 December, 2012

11981940 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.