

Irisin (42-112) (Human, Rat, Mouse, Canine)

Catalog#	Standard size	Price	Order
067-17	100 µg	\$350.00	<input type="button" value="ORDER"/>

Product description

More information

[Irisin - Antibody Stainings, Sequences & References](#)

[Irisin Synthesis and Its Applications](#)

Catalog #	067-17
Standard Size	100 µg
Sequence	Asp - Val - Arg - Met - Leu - Arg - Phe - Ile - Gln - Glu - Val - Asn - Thr - Thr - Thr - Arg - Ser - Cys - Ala - Leu - Trp - Asp - Leu - Glu - Glu - Asp - Thr - Glu - Tyr - Ile - Val - His - Val - Gln - Ala - Ile - Ser - Ile - Gln - Gly - Gln - Ser - Pro - Ala - Ser - Glu - Pro - Val - Leu - Phe - Lys - Thr - Pro - Arg - Glu - Ala - Glu - Lys - Met - Ala - Ser - Lys - Asn - Lys - Asp - Glu - Val - Thr - Met - Lys - Glu
Species	Human, Mouse, Rat
M.W	8170.25
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Insoluble in water, Dissolve peptide in 20% acetonitrile or small amount of DMSO and then dilute with water or buffer of pH>8.0
Storage Condition	Store peptide in DRY form at 0-5°C For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions. Peptide contains free cysteine and tends to be dimerized.
Appearance	White powder
Content	100 µg NET peptide
HPLC Analysis	Purity ≥ 90%

Thursday 20 December, 2012

11981942 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.