

FNDC5, Isoform 4 (Human, Rat, Mouse)

| Catalog# | Standard size | Price | Order |
|----------|---------------|----------|--------------------------------------|
| 067-18 | 10 µg | \$350.00 | <input type="button" value="ORDER"/> |

Product description

More information

[Irisin - Antibody Stainings, Sequences & References](#)

| | |
|-------------------|--|
| Catalog # | 067-18 |
| Standard Size | 10 µg |
| Sequence | MLRFI QEVNT TTRSC ALWDL EEDTE YIVHV QAISI QGQSP ASEPV LFKTP REAEK MASKN KDEVT MKEMG RNQQL RTGEV LIIVV VLFMW AGVIA LFCRQ YDIIK DNEPN NNKEK TKSAS ETSTP EHQQG GLLRS KI |
| Species | Human, Rat, Mouse |
| Storage Condition | Store in dry form at -80°C. For best results and reproducibility, rehydrate protein just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions. |
| Appearance | White powder |
| Content | 10 µg NET peptide |

Thursday 20 December, 2012

11981943 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.