

Quality Control Record

Catalog No.: B-007-25

Lot No.: 430358

Product: Biotinyl-ELA-32 (Human)

Sequence: Biotinyl--Gln-Arg-Pro-Val-Asn-Leu-Thr-Met-Arg-Arg-Lys-Leu-Arg-Lys-His-Asn-Cys-Leu-Gln-Arg-Arg-Cys-Met-Pro-Leu-His-Ser-Arg-Val-Pro-Phe-Pro
Disulfide bridge between Cys¹⁷ and Cys²²

Molecular Weight:
4194.15

Mass Spec Analysis:
Exhibits correct M.W.

HPLC Analysis:
Purity \geq 90%

Appearance: White powder

Peptide Content:
This vial contains 20 μ g of NET peptide.

Solubility and Reconstitution:
Soluble in water.

Storage Conditions:
Store peptide in DRY form at 0-5°C.

Notes: For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions.

Quality Control Manager: _____ **Date:** _____
Galina Citowsky