

Apelin-12 (Human, Rat, Mouse, Bovine) - Chemiluminescent EIA Kit

Catalog#	Standard size	Price	Order
CEK-057-23	1 kit	\$460.00	ORDER

Product description

More information

[Apelin - Research Abstracts, Antibody Stainings, Kits Standard Curve, Sequences, etc.](#)


Normal Plasma Peptide/Protein Levels

Protocol	Assay Protocol (PDF)
Sample Preparation	Blood collection Plasma extraction Tissue extraction CSF extraction
Catalog #	CEK-057-23
Standard Size	96 Wells
Sensitivity	Minimum Detectable Concentration = 13.7 pg/ml
Range	0 - 10000 pg/ml
Linear Range	13.7-151 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Apelin-12 (Human)	100
	Apelin-13 (Human)	100
	Apelin-13 (Rat)	100
	Apelin-36 (Human)	100
	Apelin-36 (Rat)	100
	Alpha-ANP (Human)	0
	ANP (25-56) (N-ANP)	0
	ADM-52 (Human)	0
	BNP-32 (Human)	0
	CNP-22	0
	Ghrelin (Human)	0
	ET-1 (Human)	0
	Bradykinin	0

Sequence Arg - Pro - Arg - Leu - Ser - His - Lys - Gly - Pro - Met - Pro - Phe

Standard
Curve Image



Thursday 20 December, 2012

11983257 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.