

Chemerin / TIG2 (145-157), Tyr145, Phe149, (Human) - EIA Kit

Catalog#	Standard size	Price	Order
EK-002-49	1 kit	\$450.00	<input type="button" value="ORDER"/>

Product description

References

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-002-49

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.27 ng/ml

Range 0-100 ng/ml

Linear Range 0.27-5.4 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Chemerin / TIG-2 (145-157)	100
	Chemerin / TIG-2 (135-154)	0
	Chemerin / TIG-2 (98-137)	0
	Chemerin / TIG-2 (72-93)	0
	Chemerin / TIG-2 (29-67)	0
	Alpha-MSH	0
	Ghrelin (Human)	0
	NPY (Human)	0
	Orexin A	0
	Secretin (Human)	0
	VIP	0

Sequence Tyr - His - Ser - Phe - Phe - Phe - Pro - Gly - Gln - Phe - Ala - Phe - Ser

Standard
Curve Image

