

# Leptin (Mouse) - EIA Kit - Extraction Free

[Leptin](#)

[Leptin Sequence Comparison](#)

**Protocol**            [Please Inquire \(email\)](#) - **Extraction Not Required**

**Sample Preparation**    [Blood collection](#)  
[Plasma extraction](#)  
[Tissue extraction](#)  
[CSF extraction](#)

**Catalog #**            EK-003-13B

**Standard Size**        96 wells

**Sensitivity**            Minimum Detectable Concentration = 0.63 ng/ml

**Range**                0.03-300 ng/ml

**Linear Range**        0.63-13ng/ml

<b>Cross Reactivity</b>	<b>Peptide</b>	<b>% Cross-reactivity</b>
	Leptin (Mouse)	100
	Leptin (Rat)	0
	Leptin (Human)	0
	TNF-alpha (Human)	0
	IL-6 (Human)	0
	IGFBP-3 (Human)	0
	IGFBP-5 (Human)	0
	IGFBP-6 (Human)	0

Tuesday 23 December, 2008

2473065 requests since Wednesday 05 April, 2006

Copyright © 2008 PHOENIX PHARMACEUTICALS, INC.