

Leptin (Rat, Mouse) - ELISA Kit

Catalog#	Standard size	Price	Order
EK-003-17	1 kit	\$495.00	<input type="button" value="ORDER"/>

- Product description**
- More information
- References

[Leptin](#)

[Leptin Sequence Comparison](#)

Protocol	Assay Protocol (PDF)	
Sample Preparation	Blood collection Tissue extraction CSF extraction	
Catalog #	EK-003-17	
Standard Size	96 wells	
Sensitivity	Minimum Detectable Concentration = 62.5 pg/ml	
Linear Range	62.5-4000 pg/ml	
Cross Reactivity	Peptide	% Cross-reactivity
	Leptin (Mouse)	100
	Leptin(Rat)	>90
	Leptin (Human)	0
	Insulin (Mouse)	0
	Somatostatin (Human)	0
	C-Peptide (Rat)	0
	PPP (Rat)	0
	Glucagon (1-37)	0
	Glucagon (1-29)	0
	IGF-I	0
	IGF-II	0

Standard Curve Image

