

# LANP / ANP (26-55), prepro (Human) - EIA Kit

Catalog#	Standard size	Price	Order
EK-005-19	1 kit	\$420.00	<input type="button" value="ORDER"/>

## Product description

## More information

ANP (Atrial Natriuretic Peptide)

Normal Plasma Peptide/Protein Levels

**Protocol** [Assay Protocol \(PDF\)](#)

**Sample Preparation** [Blood collection](#)  
[Plasma extraction](#)  
[Tissue extraction](#)  
[CSF extraction](#)

**Catalog #** EK-005-19

**Standard Size** 96 wells


**Sensitivity** Minimum Detectable Concentration = 0.32ng/ml

**Range** 0-100ng/ml

**Linear Range** 0.32-7.6ng/ml

<b>Cross Reactivity</b>	<b>Peptide</b>	<b>%Cross-reactivity</b>
	Prepro-ANP (26-55) (Human)	100
	$\alpha$ -ANP (1-28) (Human, Canine, Porcine, Bovine)	0
	ANP (Rat)	0
	Prepro-ANP (56-92) (Human)	0
	BNP-32 (Human)	0
	CNP-22 (Human, Porcine, Rat)	0
	Insulin (Human)	0
	IGF-1 (Human)	0

**Sequence** Asn - Pro - Met - Tyr - Asn - Ala - Val - Ser - Asn - Ala - Asp - Leu - Met - Asp - Phe - Lys - Asn - Leu - Leu - Asp - His - Leu - Glu - Glu - Lys - Met - Pro - Leu - Glu - Asp



○ )