

ANP (Rat, Mouse) - EIA Kit

Catalog#	Standard size	Price	Order
EK-005-24	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

References

ANP (Atrial Natriuretic Peptide)

Normal Plasma Peptide/Protein Levels

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-005-24

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.2 ng/ml

Range 0-100 ng/ml

Linear Range 0.2-5.5ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	ANP (1-28) (Rat, Mouse)	100
	α -ANP (1-28) (Human, Canine, Porcine, Bovine)	100
	α -ANP (8-33) (Rat)	100
	β -ANP (1-28) (Human)	100
	Urodilatin	100
	BNP-26 (Porcine)	<0.1
	Endothelin-1 (Human, Rat, Porcine, Mouse, Canine, Bovine)	0
	BNP-32 (Human)	0
	BNP-45 (Rat)	0
	CNP-22 (Human, Porcine, Rat)	0

Sequence Ser - Leu - Arg - Arg - Ser - Ser - Cys - Phe - Gly - Gly - Arg - Ile - Asp - Arg - Ile - Gly - Ala - Gln - Ser - Gly - Leu - Gly - Cys - Asn - Ser - Phe - Arg - Tyr

Disulfide Bridge: Cys7 - Cys23

