

Bradykinin (Human, Rat, Mouse) - EIA Kit

Catalog#	Standard size	Price	Order
EK-009-01	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

References

Bradykinin and Related Peptides

Normal Plasma Peptide/Protein Levels

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-009-01

Standard Size 96 wells


Sensitivity Minimum Detectable Concentration = 0.09 ng/ml

Range 0-100 ng/ml

Linear Range 0.09-1.13 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Bradykinin	100
	[Lys]-Bradykinin	100
	Biotinyl-Bradikinin	100
	[Tyr ⁰]-Bradykinin	100
	[Des-Arg ⁰]-Bradykinin	0
	Bradykinin Potentiator B	0
	Angiotensin II (Human)	0
	Met-Enkephalin	0
	Substance P	0
	Endothelin-1 (Human, Porcine, Rat, Mouse, Canine, Bovine)	0
	TRH	0

Sequence Arg - Pro - Pro - Gly - Phe - Ser - Pro - Phe - Arg



o