

Adrenomedullin N-20, Pro / PAMP-20 / Prodepin (Human) - EIA Kit

Catalog#	Standard size	Price	Order
EK-010-05	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

[Adrenomedullin](#)

[Normal Plasma Peptide/Protein Levels](#)

[Adrenomedullin-5 Porcine](#)

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-010-05

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.08 ng/ml

Range 0-100 ng/ml

Linear Range 0.08-0.6 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Pro-Adrenomedullin N-20 (Human)	100
	Pro-Adrenomedullin N-20 (Rat)	0
	Adrenomedullin (1-52) (Human)	0
	Adrenomedullin (1-50) (Rat)	0
	Amylin (Human)	0
	Endothelin-1 (Human, Porcine, Rat, Mouse, Canine, Bovine)	0

Sequence Ala - Arg - Leu - Asp - Val - Ala - Ser - Glu - Phe - Arg - Lys - Lys - Trp - Asn - Lys - Trp - Ala - Leu - Ser - Arg - NH₂

