

Adrenomedullin (1-50) (Rat) - EIA Kit

Catalog#	Standard size	Price	Order
EK-010-08	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

References

Adrenomedullin

Normal Plasma Peptide/Protein Levels

Adrenomedullin-5 Porcine

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-010-08

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.1 ng/ml

Range 0-100 ng/ml

Linear Range 0.1-1.79 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Adrenomedullin (Rat)	100
	Adrenomedullin (Human)	100
	Adrenomedullin (1-50) (Mouse)	100
	Pro-ADM N20 (Human)	0
	Pro-ADM N20 (Rat)	0
	Amylin (Human)	0
	Endothelin-1 (ET-1) (Human, Bovine, Canine, Mouse, Porcine, Rat)	0

Sequence Tyr - Arg - Gln - Ser - Met - Asn - Gln - Gly - Ser - Arg - Ser - Thr - Gly - Cys - Arg - Phe - Gly - Thr - Cys - Thr - Met - Gln - Lys - Leu - Ala - His - Gln - Ile - Tyr - Gln - Phe - Thr - Asp - Lys - Asp - Lys - Asp - Gly - Met - Ala - Pro - Arg - Asn - Lys - Ile - Ser - Pro - Gln - Gly - Tyr - NH₂

Disulfide Bridge: Cys14 - Cys19

