

# BNP-32 (Human) - EIA Kit

Catalog#	Standard size	Price	Order
EK-011-03	1 kit	\$420.00	<input type="button" value="ORDER"/>

## Product description

## More information

## References

Brain Natriuretic Peptides (BNP)

Normal Plasma Peptide/Protein Levels

**Protocol** [Assay Protocol \(PDF\)](#)

**Sample Preparation** [Blood collection](#)  
[Plasma extraction](#)  
[Tissue extraction](#)  
[CSF extraction](#)

**Catalog #** EK-011-03

**Standard Size** 96 wells

**Sensitivity** Minimum Detectable Concentration = 0.26ng/ml



**Range** 0-100 ng/ml

**Linear Range** 0.26-3.33 ng/ml

<b>Cross Reactivity</b>	<b>Peptide</b>	<b>%Cross-reactivity</b>
	BNP-32 (Human)	100
	BNP-32 (Rat)	<0.04
	BNP-45 (Rat)	0
	BNP-26 (Porcine)	0
	$\alpha$ -ANP (1-28) (Human, Canine)	0
	Endothelin-1 (Human, Rat, Porcine, Mouse, Canine, Bovine)	0
	Angiotensin II (Human)	0
	VIP (Human, Rat, Porcine)	0

**Sequence** Ser - Pro - Lys - Met - Val - Gln - Gly - Ser - Gly - Cys - Phe - Gly - Arg - Lys - Met - Asp - Arg - Ile - Ser - Ser - Ser - Ser - Gly - Leu - Gly - Cys - Lys - Val - Leu - Arg - Arg - His

**Disulfide Bridge:** Cys10 - Cys26

Standard  
Curve Image

