

# BNP-26 (Porcine, Ovine) - EIA Kit

| Catalog#  | Standard size | Price    | Order                                |
|-----------|---------------|----------|--------------------------------------|
| EK-011-08 | 1 kit         | \$420.00 | <input type="button" value="ORDER"/> |

## Product description

## More information

Brain Natriuretic Peptides (BNP)

Normal Plasma Peptide/Protein Levels

**Protocol** [Assay Protocol \(PDF\)](#)

**Sample Preparation** [Blood collection](#)  
[Plasma extraction](#)  
[Tissue extraction](#)  
[CSF extraction](#)

**Catalog #** EK-011-08

**Standard Size** 96 wells


**Sensitivity** Minimum Detectable Concentration = 0.27 ng/ml

**Range** 0-100 ng/ml

**Linear Range** 0.27-3.44 ng/ml

| Cross Reactivity | Peptide                          | %Cross-reactivity |
|------------------|----------------------------------|-------------------|
|                  | BNP-26 (Porcine, Ovine)          | 100               |
|                  | BNP-32 (Porcine)                 | 100               |
|                  | BNP-32 (Human)                   | 0                 |
|                  | BNP-32 (Rat)                     | 0                 |
|                  | BNP-45 (Rat)                     | 0                 |
|                  | alpha-ANP (1-28) (Human, Canine) | 0                 |

Standard Curve Image



Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.