

Augurin, Prepro (71-107) (Human) - EIA Kit

Catalog#	Standard size	Price	Order
EK-012-22	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

[Augurin - Research Abstracts, Antibody Stainings, Kits Standard Curves, Sequences, etc.](#)

[Isolation and Identification of Augurin Peptide from the Bovine Hypothalamus](#)

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-012-22

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.16ng/ml

Range 0-100 ng/ml

Linear Range 0.16-1.46 ng/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Prepro-Augurin (71-107) (H)	100
	Prepro-Augurin (108-132) (H)	0
	Prepro-Augurin (133-148) (H)	0

Sequence pGlu - Leu - Trp - Asp - Arg - Thr - Arg - Pro - Glu - Val - Gln - Gln - Trp - Tyr - Gln - Gln - Phe - Leu - Tyr - Met - Gly - Phe - Asp - Glu - Ala - Lys - Phe - Glu - Asp - Asp - Ile - Thr - Tyr - Trp - Leu - Asn - Arg

Standard Curve Image



