

# Augurin, Prepro (108-132) (Human) - EIA Kit

Catalog#	Standard size	Price	Order
EK-012-23	1 kit	\$420.00	<input type="button" value="ORDER"/>

## Product description

## More information

[Augurin - Research Abstracts, Antibody Stainings, Kits Standard Curves, Sequences, etc.](#)

**Protocol** [Assay Protocol \(PDF\)](#)

**Sample Preparation** [Blood collection](#)  
[Plasma extraction](#)  
[Tissue extraction](#)  
[CSF extraction](#)

**Catalog #** EK-012-23

**Standard Size** 96 wells

**Sensitivity** Minimum Detectable Concentration = 0.25ng/ml

**Range** 0-100 ng/ml

**Linear Range** 0.25-3.31 ng/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Prepro-Augurin (108-132) (H)	100
	Prepro-Augurin (71-107) (H)	0
	Prepro-Augurin (133-148) (H)	0

**Sequence** Asp - Arg - Asn - Gly - His - Glu - Tyr - Tyr - Gly - Asp - Tyr - Tyr - Gln - Arg - His - Tyr - Asp - Glu - Asp - Ser - Ala - Ile - Gly - Pro - Arg

**Standard Curve Image**

