

DNP (Green Mamba) - EIA Kit

Catalog#	Standard size	Price	Order
EK-013-01	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

DNP (Dendroaspis Natriuretic Peptide)

Normal Plasma Peptide/Protein Levels

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-013-01

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.42 ng/ml

Range 0-100 ng/ml

Linear Range 0.42-5.9 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	DNP	100
	α -ANP (1-28) (Human, Canine)	0
	BNP-32 (Human)	0
	CNP-22 (Human, Porcine, Rat)	0
	Prepro-ANP (26-55) (Human)	0

Sequence Glu - Val - Lys - Tyr - Asp - Pro - Cys - Phe - Gly - His - Lys - Ile - Asp - Arg - Ile - Asn - His - Val - Ser - Asn - Leu - Gly - Cys - Pro - Ser - Leu - Arg - Asp - Pro - Arg - Pro - Asn - Ala - Pro - Ser - Thr - Ser - Ala

Disulfide Bridge: Cys7 - Cys23

Standard
Curve Image

