

Amylin (Human) - EIA Kit

Catalog#	Standard size	Price	Order
EK-017-03	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

References

Amylin

Normal Plasma Peptide/Protein Levels

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-017-03

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.58ng/ml

Range 0-100 ng/ml

Linear Range 0.58-5.0ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Amylin (Human)	100
	Amylin-Amide (Human)	100
	Amylin (Rat)	0
	Amylin (Feline)	0
	CGRP (Rat)	0
	CGRP (Human)	1
	Calcitonin (Human)	0
	NPY (Porcine)	0
	Pancreastatin (Porcine)	0
	VIP (Human, Rat, Porcine)	0
	Insulin (Human)	0
	Glucagon (Human)	0
	Somatostatin-28	0
	Pancreatic Polypeptide (Human)	0
	Pro-Insulin (Human)	0

Sequence Lys - Cys - Asn - Thr - Ala - Thr - Cys - Ala - Thr - Gln - Arg - Leu - Ala - Asn - Phe - Leu - Val - His - Ser - Ser - Asn - Asn - Phe - Gly - Ala - Ile - Leu - Ser - Ser - Thr - Asn - Val - Gly - Ser - Asn - Thr - Tyr

Disulfide Bridge: Cys2 - Cys7

