

Clusterin (Human) ELISA Kit

Catalog#	Standard size	Price	Order
EK-018-35	1 kit	\$420.00	<input type="button" value="ORDER"/>
Product description			
More information			
References			

Predict Alzheimer's Disease up to 10 years in advance

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-018-35

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 1.56 ng/ml

Linear Range 1.56-100 ng/ml

Cross Reactivity	Peptide	% Cross-Reactivity
	Clusterin (Human)	100

Standard Curve Image

