

Urocortin (Human) - EIA Kit

| | | | |
|-----------|---------------|----------|--------------------------------------|
| Catalog# | Standard size | Price | Order |
| EK-019-14 | 1 kit | \$420.00 | <input type="button" value="ORDER"/> |

Product description

More information

References

[Stresscopin \(SCP\) & Stresscopin-related Peptide \(SRP\)](#)

[Urocortin Precursor Schematic \(Human\)](#)

[Urocortin and Related Peptides Sequence Comparison](#)

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-019-14

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration: 0.3ng/ml

Range 0-100ng/ml

Linear Range 0.3-3.0 ng/ml

| Cross Reactivity | Peptide | %Cross-reactivity |
|-------------------------|-------------------------------------|--------------------------|
| | Urocortin (Human) | 100 |
| | Urocortin (Rat) | 100 |
| | Urocortin II (Human) | 0 |
| | Urocortin II (Mouse) | 0 |
| | Urocortin III (Human) | 0 |
| | Urocortin III (Mouse) | 0 |
| | CRF (Human, Rat) | 0 |
| | Cortistatin14 (Human) | 0 |
| | MCH (Human, Rat, Mouse) | 0 |
| | LH-RH (Human, Porcine, Rat) | 0 |
| | NPY (Rat, Rat) | 0 |
| | Somatostatin-28 (Human, Rat, Mouse) | 0 |

Sequence Asp - Asn - Pro - Ser - Leu - Ser - Ile - Asp - Leu - Thr - Phe - His - Leu - Leu - Arg - Thr - Leu - Leu - Glu - Leu - Ala - Arg - Thr - Gln - Ser - Gln - Arg - Glu - Arg - Ala - Glu - Gln - Asn - Arg - Ile - Ile - Phe - Asp - Ser - Val - NH₂

Standard
Curve Image

