

Urocortin (Rat, Mouse) - EIA Kit

Catalog#	Standard size	Price	Order
EK-019-15	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

References

[Stresscopin \(SCP\) & Stresscopin-related Peptide \(SRP\)](#)

[Urocortin and Related Peptides Sequence Comparison](#)

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-019-15

Standard Size 96 wells


Sensitivity Minimum Detectable Concentration = 0.38 ng/ml

Range 0-100 ng/ml

Linear Range 0.38-3.47 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Urocortin (Rat)	100
	Urocortin (Human)	100
	Cortistatin-14 (Rat)	1.94
	CRF (Human, Rat)	0
	MCH (Human, Rat, Mouse)	0
	LH-RH (Human, Porcine, Rat)	0
	NPY (Human, Rat)	0
	Somatostatin	0

Sequence Asp - Asp - Pro - Pro - Leu - Ser - Ile - Asp - Leu - Thr - Phe - His - Leu - Leu - Arg - Thr - Leu - Leu - Glu - Leu - Ala - Arg - Thr - Gln - Ser - Gln - Arg - Glu - Arg - Ala - Glu - Gln - Asn - Arg - Ile - Ile - Phe - Asp - Ser - Val - NH₂



o