

Nociceptin / Orphanin FQ (Human, Rat, Mouse, Ox) - EIA Kit

Catalog#	Standard size	Price	Order
EK-021-55	1 kit	\$420.00	<input type="button" value="ORDER"/>

- Product description**
- More information
- References

[Orphanin FQ {Nociceptin} Prohormone Schematic \(Human\)](#)

[Nociceptin Sequence Comparison](#)

[Normal Plasma Peptide/Protein Levels](#)

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation
[Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-021-55

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.18 ng/ml

Range 0-100ng/ml

Linear Range 0.18-1.76 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Orphanin FQ (OFQ)	100
	Met-Enkephalin	0
	Dynorphin A	0
	β-Endorphin (Human)	0

Sequence Phe - Gly - Gly - Phe - Thr - Gly - Ala - Arg - Lys - Ser - Ala - Arg - Lys - Leu - Ala - Asn - Gln

Image



