

Endothelin-1 (ET-1) (Human, Rat, Mouse, Porcine, Bovine, Canine, Rabbit, Monkey) - EIA Kit

Catalog#	Standard size	Price	Order
EK-023-01	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

References

[Endothelin and Endothelin-like Peptides](#)

[Preproendothelin-1: Placenta, Lung, Umbilical Vein \(Human\)](#)

[Endothelin-1 Precursor \(ET-1\) Sequence](#)

[Preproendothelin-1 \(Porcine\)](#)

[Normal Plasma Peptide/Protein Levels](#)

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-023-01

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.13ng/ml

Range 0-25 ng/ml

Linear Range 0.13-1.72ng/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Endothelin-1 (Human, Porcine, Rat, Mouse, Canine, Bovine, Rabbit, Monkey)	100
	Endothelin-2 (Human, Porcine, Rat, Mouse, Canine, Bovine)	100
	Endothelin-3 (Human, Porcine, Rat, Mouse, Canine, Bovine)	100
	Big Endothelin-1 (Human)	16
	Big Endothelin-2 (Human)	7
	Angiotensin II (Human)	0
	[Arg ⁸]-Vasopressin	0

Sequence Cys - Ser - Cys - Ser - Ser - Leu - Met - Asp - Lys - Glu - Cys - Val - Tyr - Phe - Cys - His - Leu - Asp - Ile - Ile - Trp

Disulfide Bridge: Cys1 - Cys15 and Cys3 - Cys11

