

Endothelin-3 (ET-3) (Human, Rat, Porcine) - EIA Kit

Catalog#	Standard size	Price	Order
EK-023-17	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

[Endothelin and Endothelin-like Peptides](#)

[Preproendothelin-3 \(Human\)](#)

[Normal Plasma Peptide/Protein Levels](#)

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-023-17

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.33 ng/ml

Range 0-25 ng/ml

Linear Range 0.33-10.2 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Endothelin-3 (Human, Rat, Porcine, Rabbit)	100
	Endothelin-1 (Human, Rat, Porcine)	0.1
	Endothelin-2 (Human, Canine)	0.1
	Big Endothelin-1 (Porcine)	0.2
	Big Endothelin-1 (Human)	0.05
	VIC (Endothelin-Beta) (Mouse)	0.02
	BNP-26 (Porcine)	<0.01
	Sarafotoxin S6b	0.001
	α -ANP (1-28) (Human, Canine)	0
	Angiotensin II (Human)	0
	[Arg ⁸]-Vasopressin	0
	CGRP (Human)	0

Sequence Cys - Thr - Cys - Phe - Thr - Tyr - Lys - Asp - Lys - Glu - Cys - Val - Tyr - Tyr - Cys - His - Leu - Asp - Ile - Ile - Trp

Disulfide Bridge: Cys1 - Cys15 and Cys3 - Cys11

Standard
Curve Image

