

# Gastrin I (Human) - EIA Kit

Catalog#	Standard size	Price	Order
EK-027-04	1 kit	\$420.00	<input type="button" value="ORDER"/>

## Product description

## More information

## References

[Gastrin Precursor Schematic \(Human\)](#)

[Normal Plasma Peptide/Protein Levels](#)

**Protocol** [Assay Protocol \(PDF\)](#)

**Sample Preparation** [Blood collection](#)  
[Plasma extraction](#)  
[Tissue extraction](#)  
[CSF extraction](#)

**Catalog #** EK-027-04

**Standard Size** 96 wells

**Sensitivity** Minimum Detectable Concentration = 0.09 ng/ml

**Range** 0-100 ng/ml

**Linear Range** 0.09-0.97 ng/ml

<b>Cross Reactivity</b>	<b>Peptide</b>	<b>%Cross-reactivity</b>
	Gastrin I (Human)	100
	Big Gastrin I (Human)	100
	CCK (26-33) (Sulfated)	3
	CCK (26-33) (Non-Sulfated)	8
	Gastrin I (Chicken)	0.5
	Pancreatic Polypeptide (Human)	0
	GRP (Human)	0
	GIP (Porcine)	0
	Amylin Amide (Human)	0

**Sequence** pGlu - Gly - Pro - Trp - Leu - Glu - Glu - Glu - Glu - Glu - Ala - Tyr - Gly - Trp - Met - Asp - Phe - NH<sub>2</sub>

