

GLP-2 (1-34) (Human) - EIA Kit

Catalog#	Standard size	Price	Order
EK-028-14	1 kit	\$395.00	<input type="button" value="ORDER"/>

Product description

More information

GLP-1 & GLP-2

Normal Plasma Peptide/Protein Levels

Neuron Function Modulation

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-028-14

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.16 ng/ml

Range 0-100 ng/ml

Linear Range 0.16-1.88 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	GLP-2 (1-34) (Human)	100
	GLP-2 (1-34) (Rat)	100
	GLP-1 (7-36) (Human)	6
	GLP-1 / Pro-Glucagon (92-128) (Human, Rat)	5
	GLP-1 (7-37) (Human)	0
	Glucagon	0
	GIP	0
	VIP	0
	Insulin	0

Standard Curve Image

