

# Resistin, Recombinant (Mouse) - ELISA Kit

Catalog#	Standard size	Price	Order
EK-028-35	1 kit	\$420.00	<input type="button" value="ORDER"/>

## Product description

## More information

## References

[Resistin](#)

[Resistin Sequence Comparison](#)

[Normal Plasma Peptide/Protein Levels](#)

**Protocol**      [Assay Protocol \(PDF\) \\*](#)

\* - For Reference Only. Dilution conditions within this protocol may differ from each lot number. Please Contact us for an updated version of the protocol

**Sample Preparation**      [Blood collection](#)  
[Tissue extraction](#)  
[CSF extraction](#)

**Catalog #**      EK-028-35

**Standard Size**      96 wells

**Sensitivity**      Minimum Detectable Concentration = 0.016 ng/ml

**Linear Range**      0.016-1 ng/ml

<b>Cross Reactivity</b>	<b>Peptide</b>	<b>%Cross-reactivity</b>
	Resistin (Mouse)	100
	Resistin (Human)	0
	Adiponectin (Mouse)	0
	Adiponectin (Human)	0
	Leptin (Mouse)	0
	Leptin (Human)	0
	Insulin (Human)	0
	GLP1 (Human)	0
	Amylin (Human)	0
RELM- $\beta$ (Mouse)	0	
RELM- $\beta$ (Human)	0	

**Standard Curve Image**

