

Growth Hormone Releasing Factor (GHRF) (Human) - EIA Kit

| Catalog# | Standard size | Price | Order |
|-----------|---------------|----------|-----------------------|
| EK-031-02 | 1 kit | \$420.00 | ORDER |

Product description

More information

[Ghrelin Receptor Antagonist \[D-lys3\] / GHRP-6](#)

[Growth Hormone-Releasing Factor \(GHRF\) Prohormone Schematic](#)

[Normal Plasma Peptide/Protein Levels](#)

[JI-38 - Potent GHRH agonist](#)

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-031-02

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.13ng/ml

Range 0-100ng/ml

Linear Range 0.13-1.82ng/ml

| Cross Reactivity | Peptide | %Cross-reactivity |
|------------------|---|-------------------|
| | Growth Hormone Releasing Factor (Human) | 100 |
| | Growth Hormone Releasing Factor (Porcine) | 100 |
| | Growth Hormone Releasing Factor (Bovine) | 100 |
| | Growth Hormone Releasing Factor (Ovine) | 100 |
| | Growth Hormone Releasing Factor (Rat) | 0 |
| | Growth Hormone Releasing Factor (1-40)(Human) | 0 |
| | Growth Hormone Releasing Factor (1-37)(Human) | 0 |
| | VIP (Human, Porcine, Rat) | 0 |
| | Growth Hormone (Human) | 0 |
| | Secretin (Human) | 0 |

Sequence Tyr - Ala - Asp - Ala - Ile - Phe - Thr - Asn - Ser - Tyr - Arg - Lys - Val - Leu - Gly - Gln - Leu - Ser - Ala - Arg - Lys - Leu - Leu - Gln - Asp - Ile - Met - Ser - Arg - Gln - Gln - Gly - Glu - Ser - Asn - Gln - Glu - Arg - Gly - Ala - Arg - Ala - Arg - Leu - NH₂

Standard
Curve Image

