

IGF Binding Protein-1 (IGFBP-1), Recombinant (Human) - ELISA Kit

Catalog#	Standard size	Price	Order
EK-033-19	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

MGF (mechano-growth factor)

IGF-I (Insulin-Like Growth Factor I)

Protocol [Assay Protocol \(PDF\)](#) *

* - For Reference Only. Dilution conditions within this protocol may differ from each lot number. Please Contact us for an updated version of the protocol

Sample Preparation [Blood collection](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-033-19

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.15ng/ml

Range 0.15-10ng/ml

Cross Reactivity	Peptide	% Cross-reactivity
	IGFBP-I (Human)	100
	IGFBP-II (Human)	0
	IGFBP-III (Human)	0
	IGF-I (Human)	0
	IGF-II (Human)	0
	IGF-I (Mouse)	0
	IGF-II (Mouse)	0
	Insulin (Human)	0

Standard Curve Image

