

IGF Binding Protein-3 (IGFBP-3), Recombinant (Human) - ELISA Kit

//-->

Product description

More information

MGF (mechano-growth factor)

IGF-I (Insulin-Like Growth Factor I)

Protocol	Please Inquire (email)	
Sample Preparation	Blood collection Tissue extraction CSF extraction	
Catalog #	EK-033-20	
Standard Size	96 wells	
Sensitivity	Minimum Detectable Concentration = 0.78ng/ml	
Range	0.78-50 ng/ml	
Cross Reactivity	Peptide	% Cross-reactivity
	IGFBP-3 (Human)	100
	IGFBP-1 (Human)	0
	IGFBP-2 (Human)	0
	IGF-I (Human)	0
	IGF-II (Human)	0
	IGF-I (Mouse)	0
	IGF-II (Mouse)	0
	Insulin (Human)	0

Wednesday 30 March, 2011

6940164 requests since Wednesday 05 April, 2006

Copyright © 2011 PHOENIX PHARMACEUTICALS, INC.