

Brain-derived Neurotrophic Factor (BDNF) (Human) - ELISA Kit

Catalog#	Standard size	Price	Order
EK-033-22	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

Brain Derived Neurotrophic Factor (BDNF)

Protocol [Please Inquire \(email\)](#)

Sample Preparation [Blood collection](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-033-22

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 7.8pg/ml

Linear Range 7.8-500 pg/ml

Cross Reactivity

Peptide	%Cross-reactivity
BDNF	100
NGF	<3
NT-3	<3
NT-4/5	<3

Standard Curve Image



