

IL-6, Recombinant (Human) - ELISA Kit

Catalog#	Standard size	Price	Order
EK-033-31	1 kit	\$420.00	<input type="button" value="ORDER"/>
Product description	More information	References	

IL-6

Protocol [Assay Protocol \(PDF\) *](#)

* - For Reference Only. Dilution conditions within this protocol may differ from each lot number. Please Contact us for more information on the protocol

Sample Preparation [Blood collection](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-033-31

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 8pg/ml

Linear Range 8-500 pg/ml

Cross Reactivity	Peptide	Cross-reactivity
	IL -6 (Human)	100
	Adiponectin (Human)	0
	Adiponectin (Mouse)	0
	Resistin (Human)	0
	Resistin (Mouse)	0
	Leptin (Human)	0
	Insulin (Human)	0
	C-Peptide (Human)	0
	GLP1 (Human)	0
	Amylin (Human)	0
	TNF- α (Human)	0
	CNTF (Human)	0

