

C-Peptide (Human) - EIA Kit

Catalog#	Standard size	Price	Order
EK-035-01	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

[C-Peptide from Proinsulin](#)

[C-Peptide from Prepro-Insulin-like Peptide](#)

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-035-01

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration: 3.4ng/ml

Range 0.1 1000 ng / ml

Linear Range 0.26-5.4ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	C-Peptide (Human)	100
	Glucagon (Human)	0
	Insulin (Human)	0
	Amylin Amide (Human)	0
	Amylin (Human)	0
	Somatostatin	0
	Pancreatic Polypeptide (Human)	0
	CGRP (Human)	0
	CGRP (Rat)	0
	Pancreastatin (Human)	0
	VIP (Human, Porcine, Rat)	0

Sequence Glu - Ala - Glu - Asp - Leu - Gln - Val - Gly - Gln - Val - Glu - Leu - Gly - Gly - Gly - Pro - Gly - Ala - Gly - Ser - Leu - Gln - Pro - Leu - Ala - Leu - Glu - Gly - Ser - Leu - Gln

Standard
Curve Image

