

# Insulin (Human) - ELISA Kit

Catalog#	Standard size	Price	Order
EK-035-06	1 kit	\$420.00	<input type="button" value="ORDER"/>

## Product description

## More information

## References

[CRP \(Human\) ELISA Kit](#)

[Insulin \(Human\) ELISA Kit](#)

[Insulin Precursor Schematic \(Human\)](#)

[Normal Plasma Peptide/Protein Levels](#)

<b>Protocol</b>	<a href="#">Please Inquire (email)</a>
<b>Sample Preparation</b>	<a href="#">Blood collection</a> <a href="#">Tissue extraction</a> <a href="#">CSF extraction</a>
<b>Catalog #</b>	EK-035-06
<b>Standard Size</b>	96 wells
<b>Sensitivity</b>	Minimum Detectable Concentration = 0.78 uIU/ml
<b>Linear Range</b>	0.78-50 uIU/ml

<b>Cross Reactivity</b>	<b>Peptide</b>	<b>%Cross-reactivity</b>
	Insulin (Human)	100
	Ghrelin (Human)	0
	C-Peptide	0
	Somatostatin	0
	Glucagon (Human)	0
	Amylin Amide (Human)	0
	Amylin (Human)	0
	Pancreatic Polypeptide (Human)	0
	Pancreastatin (Human)	0
	IGF-1 (Human)	0
	GHRF (Human)	0
	EGF (Human)	0

Standard  
Curve Image

