

Neurotensin (Human, Rat, Mouse) - EIA Kit

Catalog#	Standard size	Price	Order
EK-048-03	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

References

Neurotensin

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-048-03

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.2 ng/ml


Range 0-100ng/ml

Linear Range 0.2-3.51 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Neurotensin	100
	Neurotensin (Guinea Pig)	26
	Neuromedin N (Porcine)	0.03
	Substance P	0
	Kinetensin	0
	Bombesin	0

Sequence pGlu - Leu - Tyr - Glu - Asn - Lys - Pro - Arg - Arg - Pro - Tyr - Ile - Leu

Standard Curve Image



Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.