

NPY (Human, Rat, Mouse) - EIA Kit

Catalog#	Standard size	Price	Order
EK-049-03	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

References

Normal Plasma Peptide/Protein Levels

Neuropeptide Y (NPY)

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-049-03

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.09 ng/ml

Range 0-100ng/ml

Linear Range 0.09-1.43 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Neuropeptide Y (Human, Rat)	100
	Neuropeptide Y (Porcine)	100
	Peptide YY (Human)	0
	Peptide YY (Rat, Mouse, Porcine, Canine)	0
	Pancreatic Polypeptide (Human)	0
	VIP (Human, Rat, Porcine)	0
	Insulin (Human)	0
	Amylin (Human)	0
	Somatostatin	0

Sequence Tyr - Pro - Ser - Lys - Pro - Asp - Asn - Pro - Gly - Glu - Asp - Ala - Pro - Ala - Glu - Asp - Met - Ala - Arg - Tyr - Tyr - Ser - Ala - Leu - Arg - His - Tyr - Ile - Asn - Leu - Ile - Thr - Arg - Gln - Arg - Tyr - NH₂

