

Chromogranin A (286-301)-Amide (Human) - EIA Kit

Catalog#	Standard size	Price	Order
EK-053-07	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

References

Chromogranin A

Serpinin

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-053-07

Standard Size 96 well

Sensitivity minimum detection level = 0.08 ng/ml

Range 0-100 ng/ml

Linear Range 0.08-1.27 ng/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Pancreastatin (Human)	100
	Pancreastatin/Chromogranin A (250-301) Amide (Human)	63
	Pancreastatin (Porcine)	28
	Pancreastatin (Rat)	3.5
	Pancreatic Polypeptide (Rat)	0
	Pancreatic Polypeptide (Human)	0
	Amylin Amide (Rat)	0
	Amylin (Human)	0
	NPY (Porcine)	0
	BNP-32 (Human)	0

Sequence Glu - Glu - Glu - Glu - Glu - Met - Ala - Val - Val - Pro - Gln - Gly - Leu - Phe - Arg - Gly - NH₂

Standard
Curve Image

