

# Pancreastatin (Porcine) - EIA Kit

Catalog#	Standard size	Price	Order
EK-053-08	1 kit	\$420.00	<input type="button" value="ORDER"/>

## Product description

## More information

[Chromogranin A](#)

[Pancreastatin](#)

[Serpinin](#)

**Protocol** [Assay Protocol \(PDF\)](#)

**Sample Preparation** [Blood collection](#)  
[Plasma extraction](#)  
[Tissue extraction](#)  
[CSF extraction](#)

**Catalog #** EK-053-08

**Standard Size** 96 wells

**Sensitivity** Minimum Detectable Concentration = 0.43 ng/ml

**Range** 0-100 ng/ml

**Linear Range** 0.43-4.12 ng/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Pancreastatin (Porcine)	100
	Pancreastatin (Human)	>100
	Pancreastatin/Chromogranin A (250-301) Amide (Human)	63
	Pancreastatin (Rat)	3.5
	Pancreatic Polypeptide (Rat)	0
	Pancreatic Polypeptide (Human)	0
	Amylin Amide (Rat)	0
	Amylin (Human)	0
	NPY (Porcine)	0
	BNP-32 (Human)	0

**Sequence** Gly - Trp - Pro - Gln - Ala - Pro - Ala - Met - Asp - Gly - Ala - Gly - Lys - Thr - Gly - Ala - Glu - Glu - Ala - Gln - Pro - Pro - Glu - Gly - Lys - Gly - Ala - Arg - Glu - His - Ser - Arg - Gln - Glu - Glu - Glu - Glu - Glu - Glu - Thr - Ala - Gly - Ala - Pro - Gln - Gly - Leu - Phe - Arg - Gly - NH<sub>2</sub>

Standard  
Curve Image

