

Chromogranin A, Prepro-(342-355) (Human) - EIA Kit

Catalog#	Standard size	Price	Order
EK-053-26	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

Chromogranin A

Serpinin

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-053-26

Standard Size 96 wells

Family Chromogranin A

Sensitivity Minimum Detectable Concentration = 0.11ng/ml


Range 0-100ng/ml

Linear Range 0.11-2.29 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Prepro-Chromogranin A (342-355) (Human)	100
	Prepro-Chromogranin A (393-426) (Human)	0
	Prepro-Chromogranin A (429-457) (Human)	0
	Prepro-Chromogranin A (392-402) (Mouse)	0

Sequence Trp - Ser - Lys - Met - Asp - Gln - Leu - Ala - Lys - Glu - Leu - Thr - Ala - Glu

Standard Curve Image



Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.