

Catestatin (Human) - EIA Kit - Extraction Free\*

Catalog#	Standard size	Price	Order
EK-053-27	1 kit	\$420.00	<input type="button" value="ORDER"/>

- Product description
- More information

Chromogranin A

Catestatin

Serpinin

Protocol      [Assay Protocol \(PDF\)](#) - **No Sample Extraction Required\***

\*Serum and plasma do not require extraction. All other sample types may need extraction.

Sample Preparation      [Blood collection](#)  
[Tissue extraction](#)  
[CSF extraction](#)

Catalog #      EK-053-27

Standard Size      96 well EIA Kit

Sensitivity      0.05 ng/ml

Range      0-100 ng/ml

Linear Range      0.05-0.92 ng/ml

Standard Curve Image

