

Catestatin (Rat) - EIA Kit - Extraction Free*

| | | | |
|-----------|---------------|----------|--------------------------------------|
| Catalog# | Standard size | Price | Order |
| EK-053-29 | 1 kit | \$420.00 | <input type="button" value="ORDER"/> |

- Product description
- More information

Chromogranin A

Catestatin

Serpinin

Protocol [Assay Protocol \(PDF\)](#) - **No Sample Extraction Required***

*Serum and plasma do not require extraction. All other sample types may need extraction.

Sample Preparation [Blood collection](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-053-29

Standard Size 96 well kit

Sensitivity 0.08ng/ml

Range 0-100 ng/ml

Linear Range 0.08-2.78 ng/ml

Standard Curve Image

