

Serpinin (Human)

Catalog#	Standard size	Price	Order
EK-053-45	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

Chromogranin A

Serpinin

Sample Preparation
[Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-053-45

Standard Size 96 Wells

Sensitivity Minimum Detectable Concentration = 0.23ng/ml

Range 0.01-100ng/ml

Linear Range 0.23-2.7ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Serpinin / Prepro-Chromogranin A (429-454) (Human)	100
	Prepro-Chromogranin A (429-457) (Human)	100
	Serpinin / Prepro-Chromogranin A (435-460) (Rat, Mouse)	0
	[pGlu]-Serpinin (Human, Rat, Mouse)	0

Standard Curve Image

