

Parathyroid Hormone-related Protein (PTH-RP) (1-34) (Human, Rat, Mouse) - EIA Kit

Catalog#	Standard size	Price	Order
EK-056-04	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description

More information

References

Parathyroid Hormone and its Related Protein

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-056-04

Standard Size 96 wells


Sensitivity Minimum Detectable concentration: 0.5ng/ml

Range 0-100 ng/ml

Linear Range 0.5-4.46 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	PTH-rP (1-34) (Human)	100
	PTH-rP (1-37) (Human)	100
	PTH-rP (7-34) Amide (Human)	100
	PTH (1-34) (Human)	0
	PTH (1-34) (Rat)	0
	[Asn ⁷⁶]-PTH 184 (Human)	0
	PTH-rP (38-64) Amide (Human)	0
	PTH-rP (67-86) Amide (Human)	0
	PTH-rP (107-138) (Human)	0
	PTH-rP (109-141) (Human)	0

Sequence Ala - Val - Ser - Glu - His - Gln - Leu - Leu - His - Asp - Lys - Gly - Lys - Ser - Ile - Gln - Asp - Leu - Arg - Arg - Arg - Phe - Phe - Leu - His - His - Leu - Ile - Ala - Glu - Ile - His - Thr - Ala



o