

Apelin-36 (Human) - EIA Kit

Catalog#	Standard size	Price	Order
EK-057-15	1 kit	\$450.00	ORDER

Product description

More information

References

[Apelin - Research Abstracts, Antibody Stainings, Kits Standard Curve, Sequences, etc.](#)

Normal Plasma Peptide/Protein Levels

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # EK-057-15

Sensitivity Minimum Detectable Concentration = 0.09ng/ml

Range 0-100ng/ml

Linear Range 0.09-1.18ng/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Apelin-36 (Human)	100
	Apelin-13 (Human)	100
	Apelin-12 (Human, Rat, Mouse, Bovine)	100
	Apelin-36 (Rat, Mouse)	100
	α-ANP (Human, Canine, Ovine, Porcine)	0
	ANP (25-56) (N-ANP)	0
	ADM-52 (Human)	0
	BNP-32 (Human)	0
	CNP-22(Human, Rat, Mouse, Porcine)	0
	Ghrelin (Human)	0
	ET-1 (Human, Rat, Mouse, Canine, Porcine, Bovine)	0
	Bradykinin (Human, Rat, Mouse)	0

Sequence Leu - Val - Gln - Pro - Arg - Gly - Ser - Arg - Asn - Gly - Pro - Gly - Pro - Trp - Gln - Gly - Gly - Arg - Arg - Lys - Phe - Arg - Arg - Gln - Arg - Pro - Arg - Leu - Ser - His - Lys - Gly - Pro - Met - Pro - Phe

